

Product Composition

Composition	Amount per serving (60ml)	Amount per 100ml
Energy	76kcal 318kJ	127kcal 530kJ
Fat	0.12g	0.2g
of which saturates	0.06g	0.1g
Carbohydrates	0.33g	0.55g
of which sugars	0.33g	0.55g
Fibre	0.2g	0.33g
Protein	20g	33.3g
Minerals		
Sodium	37mg (1.61 mmol)	61.7mg (2.68 mmol)
Potassium	13mg (0.33 mmol)	21.7mg (0.56 mmol)
Phosphorus	2.6mg (0.08 mmol)	4.3mg (0.13 mmol)
Calcium	3mg (0.08 mmol)	5mg (0.13 mmol)

Request a **FREE** sample pack of any Nutrinovo product – direct to you or direct to your patient.

- Sample packs can be delivered direct to you or direct to your patient in the community.
- The majority of sample packs are despatched within 24 hours and we aim to despatch 100% within 48 hours, meaning you or your patients are able to try any ProSource product when you need it most.

**FREE
SAMPLE
SERVICE**

Request a **FREE** sample pack, visit www.nutrinovo.com
email info@nutrinovo.com or call **01249 691301**



nutrinovo
simply innovative nutrition

Nutrinovo Limited | +44 (0)1249 691301

www.nutrinovo.com @Nutrinovo

ONE CONNECTION 20g PROTEIN



20g protein | 60ml pouch | ENFit® compatible

nutrinovo
simply innovative nutrition

Nutrino are delighted to announce the launch of ProSource TF ENFit® - the world's first ENFit compatible modular protein.

Designed specifically with both the patient and the health professional in mind, NEW ProSource TF ENFit delivers 20 grams of high-quality protein directly to the patient's feeding tube via one simple, safe connection.

- ✓ Unique ENFit compatible pouch
- ✓ 20g TF protein in just 60ml volume
- ✓ Unrivalled essential amino acid profile
- ✓ Unrivalled branch chain amino acid profile
- ✓ Very low in electrolytes
- ✓ Halal certified
- ✓ Protein Digestibility Score of 100% (PDCAAS)

An unrivalled amino acid profile that guarantees you are giving your compromised patients the most effective protein, coupled with a delivery system which saves on both nursing time and costs, makes ProSource TF ENFit the ideal protein to use with your enterally fed patients.

Unrivalled Essential Amino Acid Profile

Amino Acid	Institute of Medicine Reference Pattern (mg per g protein)	ProSource TF ENFit (mg per g protein)	ProSource TF ENFit (grams per 60ml serving)
Histidine	17	18	0.36
Isoleucine	23	24	0.48
Leucine	52	55	1.1
Lysine	47	51	1.02
Methionine & Cysteine	23	24	0.48
Phenylalanine & Tyrosine	41	44	0.88
Threonine	26	28	0.56
Tryptophan	6	6	0.12
Valine	29	31	0.62

ONE STEP CONNECTION

It's never been easier to provide additional protein to your tube-fed patients – simply connect and squeeze to deliver 20 grams of ProSource TF direct into your patient's enteral feeding tube.

ProSource TF ENFit's one-step connection saves on nursing time, and omits the need to use additional syringes for administration.

DIRECT PROTEIN DELIVERY

